

BRODA[®]

Encore

Rehab Wheelchair



mobility
product award 2021

HCPCS
E1161



Back canes fold down for easy storage options



Dynamic rocking function soothes and stimulates users



Front pivot seat tilt maintains a proper position for self-propulsion



Open seating system allows for almost any seating configuration

Shown with Matrix Libra cushion, Matrix E2 back, and Matrix Loxx Head support

Patient, Facility, and Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Allows for easy self-mobility and safer patient maneuverability
- Reduces risk of injury and falls with supportive stable positioning

brodaseating.com | 1-844-552-7632

matrix[®]

SEATING SERIES



Loxx Headrest

Highly adjustable hardware with tool-less adjustment levers for easy fitting and position changes.



E2 Back

Shell - Contoured aluminum with cutouts for lifting and breathability. Available in standard (3") OR deep (6") contours.

Foam - HR (high resilient) foam enclosed in cover.



Libra Cushion



Foam - HR (high resilient) molded foam made with UltraFresh antimicrobial for odor protection. Contoured pelvic well, waterfall front edge.

MOBILITY WITH COMFORT,
FUNCTIONALITY, &
CUSTOMIZED SEATING

SPECIFICATIONS

Degrees of Tilt	Up to 33° (with Dynamic Rocking)
Weight Capacity	Up to 250 lbs. (114 kg)
Seat Widths	16", 18", 20", 22" (40cm, 45cm, 50cm, 55cm)
Seat Depths	16", 17", 18", 19", 20" (40cm, 43cm, 45cm, 48cm, 50cm)
Back Cane Heights	16", 18", 21", 24" (40cm, 45cm, 53cm, 60cm) Angle adjustable to 5°, 10°, 15°, 20°
Adjustable Floor-to-Seat Heights	14", 15", 16", 17", 18", 19" (35cm, 38cm, 40cm, 43cm, 45cm, 48cm)
Wheels	20" or 22" Front Mag (50cm or 55cm) 5" Rear Casters: Total Lock (12.5cm) 5" Rear Aluminum: Non-Locking (12.5cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame

FRAME OPTIONS



Additional accessories available: foot supports, trays, and positioning straps

brodaseating.com | 1-844-552-7632